



Video Intubation Laryngoscopy Emergency, Ambulance, ICU

التنظير الحنجري بالفيديو
للتنبيب (الطوارئ - الإسعاف - العناية المركزة)


Arabian Blue Diamond
Video Intubation Laryngoscopy



TRACHWAY®
VIDEO INTUBATION SERIES

Video Intubation Laryngoscopy

التنظير الحنجري بالفيديو للتببيب



Video Laryngoscope

Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752



**Click for
more details**

TRACHWAY®
VIDEO INTUBATION SERIES

Intubation Handel Set

Trachway® Handle Set TVI-4000



60°


45°

Adjustable angle

Thumb operation

60° upward

45° downward



Video output port



Micro SD card slot



Mini USB Charging port



TVI-4056-B Aluminum carrying case
L 508 W 270 H 100 (mm)
Weight 1.5kg

Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752



Click for
more details

TRACHWAY®
VIDEO INTUBATION SERIES

Video Malleable Stylet



► Video Light Stylet


Light Source	Surround LEDs Guiding red LED (for Ø 5.0)
Material	S.S. 304
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4102 Ø 5.0 L 345mm
- TVI-4105 Ø 5.0 L 420mm
- TVI-4106 Ø 4.0 L 300mm
- TVI-4108 Ø 4.0 L 420mm



Malleability

Bend for difficult airway case and nasal intubation



Advantage

- Patented guiding red LED
- 304 stainless steel
- Submersible for disinfection
- Seamless seal and foggy free



Quick switch connector

Click-in design for quick switch among blade, stylet, and soft style systems



Four Sizes

TVI-4102 Ø 5.0mm L: 345mm
-General ETT (5.5 above)
-Reinforce ETT

TVI-4105 Ø 5.0mm L: 420mm
-35 Fr. DLT, Preform ETT

TVI-4106 Ø 4.0mm L: 300mm
-Pediatric ETT (4.5 above)

TVI-4108 Ø 4.0mm L: 420mm
-32 Fr. DLT, Preform ETT

Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752



**Click for
more details**

TRACHWAY®
VIDEO INTUBATION SERIES

Video Flexible Stylet



► **Video Flexible Stylet**

Light Source	Surround LEDs
Material	Medical Grand PU
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4301 Ø 5.0 L 650mm
- TVI-4302 Ø 4.0 L 750mm



Flexible Tube

Not easy to break



Plug in & out connector

Fast and efficient
when implement
tube exchange



Two Sizes

TVI-4301 OD 5.0mm Lenth 750mm
TVI-4302 OD 4.0mm Lenth 750mm

Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752



Click for
more details

TRACHWAY
VIDEO INTUBATION SERIES

Video Rigid Laryngeal Blades



Video Flexible Blade

Light Source	LEDs
Material	S.S. 316
Focus Range	5-50mm
Field of view	FOV (D):83°

• TVL-4220	MILL #0	Neonatal
• TVL-4221	MILL #1	Infant
• TVL-4232	MAC #2	Pediatric
• TVL-4233	MAC #3	Adult
• TVL-4234	MAC #4	Adult



Durability

- Submersible for disinfection
- Seamless seal
- Foggy free



Advantage

- Conventional blade design
- Same skill with non direct view
- 316 stainless steel



Five Sizes

TVL-4220 Mill#0	Neonatal
-W: 11.5mm L: 75mm	
TVL-4221 Mill#1	Infant
-W: 11.5mm L: 102mm	
TVL-4232 MAC#2	Pediatric
-W: 12mm L: 100mm	
TVL-4233 MAC#3	Adult
-W: 15mm L: 130mm	
TVL-4234 MAC#4	Adult
-W: 15mm L: 135mm	



Quick switch connector

Click-in design for quick switch among blade, style, and soft style systems



Arabian Blue Diamond

Al Faisaliah Tower, Olaya District, Riyadh
P.O Box 12234, Riyadh 3124 Kingdom of Saudi Arabia
Tel: +966 (11) 490 3741
Fax: +966 (11) 490 3752



Click for
more details