



## Surgical Endoscopes All Medical Specialties

المناظير الجراحية الألمانية  
لكل التخصصات الطبية

ABD Mideast  
Surgical Endoscopes - All Specialty -Emos Germany



## Endoscopes Different Specialties

مناظير لتخصصات طبية مختلفة



### ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details




### Advantages of flexible Emos endoskopes



**Modern design**



**Made in Germany**



**Brilliant image representation**

High-quality components „made in Germany“  
 High depth of field and homogeneous illumination



**Stable guiding element**

Increasing lifetime due to innovative construction method



**Robust shaft**

Very good material resistance during the preparation process



**Focusable**

Optimal image setting through the ergonomic focusing ring



**Standard eye piece**

for camera application. Compatible with different camera TV coupler



**Universal high guide / cable connection**

Compatible with adapters from different manufacturers.

### Advantages of Emos optics



**Sapphire negative & sapphire window**

Solded sapphire negative and sapphire window  
 advantages: no glue, no leakage and no scratches



**Locking**

Compatible with sheaths from different manufacturer  
 advantages: easy combination with your existing equipment



**Autoclavable**

134°C / 273 °F / 2,3 bar / 21 min



**Material**

Complete made from stainless steel for high corrosion resistance



**3-tube-system**

Optic system is additional protected



**Laser welding**

Trough the laser welding absolutely sealed



**Image quality**

HD image display in real reproduction of the colors and tissue structures



**Image quality**

High-resolution, vignetting-free and distortion-free image quality in 4K

### ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
 Dubai, UAE PO Box 31303  
 Tel: +971 4 319 7867  
 Fax: +971 4 319 7897



**Click for  
 more details**

## ADVANTAGES OF FLEXIBLE EMOS ENDOSCOPES *Good to know*



MODERN DESIGN

MADE IN GERMANY



BRILLIANT IMAGE REPRESENTATION

High-quality components „made in Germany“  
High depth of field and homogeneous illumination



STABLE GUIDING ELEMENT

Increasing lifetime due to innovative construction method



ROBUST SHAFT

Very good material resistance during the preparation process



FOKUSABLE

Optimal image setting through the ergonomic focusing ring



STANDARD EYE PIECE for camera application

Compatible with different camera TV coupler



UNIVERSAL LIGHT GUIDE / CABLE CONNECTION

Compatible with adapters from different manufacturers

### ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details

## Bronchoscopy

مناظير القصبات الهوائية



emos  
technology



Manufacturer

Origin

Product

Arabic Name

System Photo

# Flexibles- Bronchoskop 5.2/6.0 Flexible-Bronchoscope 5.2/6.0



## Technische Daten Technical Data

| Model             | / model              | FBS-5223                                                                                         | FBS-6028                                                                                         |
|-------------------|----------------------|--------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| Einführschlauch   | / insertion tube     | Ø 5,2mm                                                                                          | Ø 6,0mm                                                                                          |
| Arbeitskanal      | / working channel    | Ø 2,3mm                                                                                          | Ø 2,8mm                                                                                          |
| Nutzlänge         | / working length     | 600mm                                                                                            | 600mm                                                                                            |
| Blickwinkel       | / field of view      | 100°                                                                                             | 100°                                                                                             |
| Abwinklung auf/ab | / angulation up/down | 180°/130°                                                                                        | 180°/130°                                                                                        |
| Zubehör           | / accessories        | Absaugventil,<br>Dichtigkeitstester,<br>Koffer<br><br>Suction valve,<br>Leak tester,<br>suitcase | Absaugventil,<br>Dichtigkeitstester,<br>Koffer<br><br>Suction valve,<br>Leak tester,<br>suitcase |
| Artikel Nr.       | / item no.           | FBS-5223C ■                                                                                      | FBS-6028C ■                                                                                      |
|                   |                      | FBS-5223B ■                                                                                      | FBS-6028B ■                                                                                      |

### ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details

Otolaryngology

مناظير الأنف والأذن والحنجرة



► Sinuskope / Otoskope /  
Laryngoskope  
Sinusscopes / Otopscopes / Laryngoscopes



► Naso-Pharyngo-  
Laryngoskop 2.8/3.2/3.4  
Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



**ABD MIDEAST GENERAL TRADING**

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details



Urology

مناظير المسالك



emos  
technology



► Semi-Flexibles-  
Ureterorenoskop  
Semi-Flexible-Ureterorenoscope



► Flexibles  
Uretero-Renoskop 8.4 Fr.  
Flexible Uretero-Renoscope 8.4 Fr.



**ABD MIDEAST GENERAL TRADING**

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details

Nephroscope

منظير الكلي



emos  
technology



# Semi-Flexibles Nephroskop

Semi-Flexible Nephroscope



## ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details



Gynecology

مناظير طب النساء



emos  
technology



## ► Zystoskope / Hysteroskope Cystoscopes / Hysteroscopes



## ► Flexibles- Zystoskop 5.5 Flexible-Cystoscope 5.5



### ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details

Laparoscopy

مناظير البطن



emos  
technology



# Laparoskop 4K

Laparoscope 4K

NEW

4K  
ENDOSCOPY



## ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details

Arthroscopy

مناظير طب العظام والمفاصل



emos  
technology



# Arthroskop Arthroscope



## ABD MIDEAST GENERAL TRADING

Level 41, Emirates Towers Sheikh Zayed Road.  
Dubai, UAE PO Box 31303  
Tel: +971 4 319 7867  
Fax: +971 4 319 7897



Click for  
more details