



الماسة العربية الزرقاء الطبية
Arabian Blue Diamond

Video Intubation
Laryngoscopy
Emergency,
Ambulance,
ICU

التنظير الحنجري
بالفيديو للتبيب
للطوارئ
للإسعاف
للغاية المركزة




A distinctive approach to life-saving..

TRACHWAY[®]
VIDEO INTUBATION SERIES


Video Intubation Laryngoscopy

التنظير الحنجري بالفيديو للتببيب



Video Laryngoscope

Trachway[®] Handle Set TVI-4000



TVI-4056-B Aluminum carrying case
L 508 W 270 H 100 (mm)
Weight 1.5kg

Video Malleable Stylet



▶ Video Light Stylet


Light Source	Surround LEDs Guiding red LED (for Ø 5.0)
Material	S.S. 304
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4102 Ø 5.0 L 345mm
- TVI-4105 Ø 5.0 L 420mm
- TVI-4106 Ø 4.0 L 300mm
- TVI-4108 Ø 4.0 L 420mm



Malleability

Bend for difficult airway case and nasal intubation



Advantage

- Patented guiding red LED
- 304 stainless steel
- Submersible for disinfection
- Seamless seal and foggy free



Quick switch connector

Click-in design for quick switch among blade, stylet, and soft style systems



Four Sizes

TVI-4102 Ø 5.0mm L: 345mm
-General ETT (5.5 above)
-Reinforce ETT

TVI-4105 Ø 5.0mm L: 420mm
-35 Fr. DLT, Preform ETT

TVI-4106 Ø 4.0mm L: 300mm
-Pediatric ETT (4.5 above)

TVI-4108 Ø 4.0mm L: 420mm
-32 Fr. DLT, Preform ETT

Video Flexible Stylet



► **Video Flexible Stylet**

Light Source	Surround LEDs
Material	Medical Grand PU
Focus Range	5-50mm
Field of view	FOV (D): 83°

- TVI-4301 Ø 5.0 L 650mm
- TVI-4302 Ø 4.0 L 750mm



Flexible Tube

Not easy to break



Plug in & out connector

Fast and efficient when implement tube exchange



Two Sizes

TVI-4301 OD 5.0mm Lenth 750mm
TVI-4302 OD 4.0mm Lenth 750mm



Video Flexible Blade

Light Source LEDs
 Material S.S. 316
 Focus Range 5-50mm
 Field of view FOV (D):83°

- TVL-4220 MILL #0 Neonatal
- TVL-4221 MILL #1 Infant
- TVL-4232 MAC #2 Pediatric
- TVL-4233 MAC #3 Adult
- TVL-4234 MAC #4 Adult



Durability

- Submersible for disinfection
- Seamless seal
- Foggy free



Quick switch connector

Click-in design for quick switch among blade, style, and soft style systems



Advantage

- Conventional blade design
- Same skill with non direct view
- 316 stainless steel



Five Sizes

- | | |
|---------------------|------------------|
| TVL-4220 Mill#0 | Neonatal |
| -W: 11.5mm L: 75mm | |
| TVL-4221 Mill#1 | Infant |
| -W: 11.5mm L: 102mm | |
| TVL-4232 MAC#2 | Pediatric |
| -W: 12mm L: 100mm | |
| TVL-4233 MAC#3 | Adult |
| -W: 15mm L: 130mm | |
| TVL-4234 MAC#4 | Adult |
| -W: 15mm L: 135mm | |