


# ▶ Sinuskope / Otoskope / Laryngoskope

Sinusscopes / Otoscopes / Laryngoscopes



*grenzenlos im Detail!* ▶▶



# Sinuskope

## Sinusscopes



Visionary Medical Technology



### Technische Daten

#### Technical Data

| Model / model              | Sinuskop Sinusscopes | Sinuskop Sinusscopes |                      |
|----------------------------|----------------------|----------------------|----------------------|
| Durchmesser / diameter     | Ø 2,8mm              | Ø 2,8mm              |                      |
| Kompatibel / compatible    | Storz®               | Storz®               |                      |
| Beschreibung / description | l=175mm / 0°         | l=175mm / 30°        |                      |
| Artikel Nr. / item no.     | SS28A                | SS28B                |                      |
| Model / model              | Sinuskop Sinusscopes | Sinuskop Sinusscopes | Sinuskop Sinusscopes |
| Durchmesser / diameter     | Ø 4mm                | Ø 4mm                | Ø 4mm                |
| Kompatibel / compatible    | Storz®               | Storz®               | Storz®               |
| Beschreibung / description | l=175mm / 0°         | l=175mm / 30°        | l=175mm / 70°        |
| Artikel Nr. / item no.     | SS40A                | SS40B                | SS40C                |
| Model / model              | Sinuskop Sinusscopes | Sinuskop Sinusscopes |                      |
| Durchmesser / diameter     | Ø 2,8mm              | Ø 2,8mm              |                      |
| Kompatibel / compatible    | Storz®               | Storz®               |                      |
| Beschreibung / description | l=110mm / 0°         | l=110mm / 30°        |                      |
| Artikel Nr. / item no.     | SS28A110             | SS28B110             |                      |
| Model / model              | Sinuskop Sinusscopes | Sinuskop Sinusscopes | Sinuskop Sinusscopes |
| Durchmesser / diameter     | Ø 4mm                | Ø 4mm                | Ø 4mm                |
| Kompatibel / compatible    | Wolf®                | Wolf®                | Wolf®                |
| Beschreibung / description | l=170mm / 0°         | l=170mm / 30°        | l=170mm / 70°        |
| Artikel Nr. / item no.     | SW40A                | SW40B                | SW40C                |





EMOS Technology GmbH  
Gewerbestr. 10  
D-88636 Illmensee

fon +49 (0) 75 58 / 93 82 78 - 0  
fax +49 (0) 75 58 / 93 82 78 - 55

E-Mail: [info@emostechnology.com](mailto:info@emostechnology.com)  
[www.emostechnology.com](http://www.emostechnology.com)

# ▶ Otokope

## Otosopes



## Technische Daten

Technical Data

| Model / model              | Otokop<br>Otoscope | Otokop<br>Otoscope |
|----------------------------|--------------------|--------------------|
| Durchmesser / diameter     | Ø 2,8mm            | Ø 2,8mm            |
| Kompatibel / compatible    | Storz®             | Storz®             |
| Beschreibung / description | l=50mm / 0°        | l=50mm / 30°       |
| Artikel Nr. / item no.     | OS28A50            | OS28B50            |
| Model / model              | Otokop<br>Otoscope | Otokop<br>Otoscope |
| Durchmesser / diameter     | Ø 4mm              | Ø 4mm              |
| Kompatibel / compatible    | Storz®             | Storz®             |
| Beschreibung / description | l=50mm / 0°        | l=50mm / 30°       |
| Artikel Nr. / item no.     | OS40A50            | OS40B50            |



# ► Laryngoskope

## Laryngoscopes



### Technische Daten

Technical Data

| Model        | / model       | Laryngoskop   | Laryngoskop   | Laryngoskop   |
|--------------|---------------|---------------|---------------|---------------|
|              |               | Laryngoscope  | Laryngoscope  | Laryngoscope  |
| Durchmesser  | / diameter    | Ø 7,0mm       | Ø 7,0mm       | Ø 4,0mm       |
| Kompatibel   | / compatible  | Universal     | Universal     | Storz®        |
| Beschreibung | / description | I=190mm / 70° | I=190mm / 90° | I=175mm / 70° |
| Artikel Nr.  | / item no.    | LU70C         | LU70D         | LS40CUD       |

# ▶▶ Naso-Pharyngo- Laryngoskop 2.8/3.2/3.4

Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



*grenzenlos im Detail!*







# Naso-Pharyngo-Laryngoskop 2.8/3.2/3.4

Naso-Pharyngo-Laryngoscope 2.8/3.2/3.4



18.000  
Pixel



## Technische Daten Technical Data

| Model                                   | / Model          | FNS-2800                                                                                      | FNS-3200                                                                                        | FNS-3400                                                                                        |
|-----------------------------------------|------------------|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| Einführschlauch                         | / insertion tube | Ø 2,8mm                                                                                       | Ø 3,2mm                                                                                         | Ø 3,4mm                                                                                         |
| Nutzlänge                               | / working length | 300mm                                                                                         | 300mm                                                                                           | 300mm                                                                                           |
| Blickwinkel                             | / field of view  | 85°                                                                                           | 85°                                                                                             | 85°                                                                                             |
| Abwinklung auf/ab<br>angulation up/down |                  | 160°/160°                                                                                     | 160°/160°                                                                                       | 160°/160°                                                                                       |
| Zubehör                                 | / accessories    | Dichtigkeitstester,<br>Transportbox<br><br>leakage tester,<br>transportbox                    | Dichtigkeitstester,<br>Transportbox<br><br>leakage tester,<br>transportbox                      | Dichtigkeitstester,<br>Transportbox<br><br>leakage tester,<br>transportbox                      |
| Artikel Nr.                             | / item no.       | FNS-2800G  | FNS-3200G  | FNS-3400G  |
|                                         |                  | FNS-2800B  | FNS-3200B  | FNS-3400B  |
|                                         |                  | FNS-2800R  | FNS-3200R  | FNS-3400R  |