


Digital otoscope

# SyncVision iO1 OTO



DESIGN  
AWARD



**Brings a reliable and easy exam**

It increases efficiency of each treatment process.  
It offers immediate communication between  
the ENT Doctors and patient.



### SyncVision iO1 OTO DELUX PACK

- One iO1 OTO digital Otoscope
- One Solid Ear Camera Probe
- One 200mm(L) Flexible Probe
- Disposable Protective Cover
- Cable Set (RCA & USB)
- AC Adapter & Wall Mount Kit



### SyncVision iO1 OTO LITE PACK

- One iO1 OTO digital Otoscope
- One Solid Ear Camera Probe
- Disposable Protective Cover
- Cable Set (RCA & USB)
- AC Adapter & Wall Mount Kit

### SyncVision iO1 OTO SPECIFICATIONS

Dimension	165 x 63 x 35 mm	Light source	High Luminance white LED
Weight	180 g	File format	JPEG
Display	2.5" full color TFT-LCD.	Interface	Video output / USB output
Probe size	Ø5mm x 50 & 120mm(L)	Memory	≥8G MicroSD
Resolution	307,200 pixels	Signal system	NTSC / PAL
Focus Type	Auto	Power	Built-in rechargeable battery (2 hours runtime)
Focus range	10~20 mm / 15~50mm		

